



[Home](#) > [Temperature](#) > Controllers

LakeShore Temperature Controllers

▶ [Selection Chart](#)



Model 370 AC Resistance Bridge



Features

- Useful down to below 20 mK
- Supports RTD sensors
- One sensor input; expandable to 16 sensor inputs
- One control loop: 1 W
- AC sinusoidal current excitation as low as 3 pico Amps
- PID temperature control
- IEEE-488 and serial interface, analog outputs, and alarm relays

[More Details](#)

[Add To Cart](#)

Model 340



Features

- Useful down to below 300 mK with appropriate RTD sensors
- Supports diode, RTD, Capacitance, and thermocouple sensors
- Two sensor inputs; expandable to ten sensor inputs
- Scalable sensor excitation current as low as 30 nA for use with NTC resistors to below 300 mK
- Sensor excitation current reversal for resistance temperature sensors
- Two control loops: 100W and 1 W
- Autotuning
- IEEE-488 and serial interface, analog outputs, digital I/O, and alarm relays

[More Details](#)

[Add To Cart](#)

Model 332





Features

- Useful down to below 1 K with appropriate RTD sensors
- Supports diode, RTD, and thermocouple sensors
- Two sensor inputs
- Scalable sensor excitation current as low as 1 uA for use with NTC resistors to below 1 K
- Sensor excitation current reversal for resistance temperature sensors
- Two control loops: 50 W and 10 W
- Autotuning
- IEEE-488 and serial interface, analog outputs and alarm relays

[More Details](#)

[Add To Cart](#)

Model 331	Model 321
	
Features <ul style="list-style-type: none"> • Useful down to 1.4 K • Supports diode, RTD, and thermocouple sensors • Two sensor inputs • Sensor excitation current as low as 10 uA for use with NTC resistors to below 1 K • Sensor excitation current reversal for resistance temperature sensors • Two control loops: 50 W and 1 W • Autotuning • IEEE-488 and serial interface, analog outputs and alarm relays 	Features <ul style="list-style-type: none"> • Useful down to 1.4 K • Supports one silicon diode, platinum RTD, or thermocouple sensor • One sensor input • Autotuning • RS-232 interface
More Details Add To Cart	More Details Add To Cart

Selection Chart							
	370	340	332	331S	331E	321	
Sensor input(s)	1 or 16	2, 4, or 10	2	2	2	1	
Control loop(s)	1	2	2	2	2	1	
Heater output(s)	1 W	100 W / 1 W	50 W / 1 W	50 W / 1 W	50 W	25 W	
RS232C interface	X	X	X	X	X	X	
IEEE-4888 interface	X	X	X	X			
Analog outout(s)	2 @ ± 10V	2 @ ± 10V	2 @ ± 10V	1 @ ± 10V		1 @ ± 10V	
Relays	X	X	X	X			
CE Mark	X	X	X	X	X		
Supported Sensors	Minimum temperature with instrument						Max*
Silicon Diodes		1.4 K	1.4 K	1.4 K	1.4 K	1.4 K	500 K
GaAlAs diodes		1.4 K	1.4 K	1.4 K	1.4 K		500 K
Cernox™	300 mK	300 mK	1 K	2 K	2 K		420 K
Germanium	< 50 mK	300 mK	1 K	2 K	2 K		100 K
Carbon-Glass™	1.4 K	1.4 K	3 K	3 K	3 K		325 K
Ruthenium Oxide	< 50 mK	300 mK	1 K	2 K	2 K		300 K
Platinum	14 K	30 K	30 K	30 K	30 K		800 K
Rhodium-iron	1.4 K	1.4 K	1.4 K	1.4 K	1.4 K	1.4 K	400 K
Capacitance		1.4 K					290 K
Thermocouples		1.4 K	1.4 K	1.4 K	1.4 K	1.4 K	1530 K

* Maximum sensor temperature with any Lake Shore temperature controller.

[Top](#)

[Home](#) | [Temperature](#) | [Magnetics](#) | [Systems](#) | [Order Now & Pricing](#) | [Contact Us](#) | [What's New](#) | [Site Map](#)

Copyright 2002 Lake Shore Cryotronics, Inc.